



## Mac OS (Classic) Installation instructions

### Welcome to CoolLocalizer

CoolLocalizer© is an image analysis software for confocal microscopy.

The program was designed to be easy to use and install. Please follow the instructions below:

### Installation

Download CoolLocalizer at [www.CoolLocalizer.com](http://www.CoolLocalizer.com).

Unpack the compressed archive using Stuffit Expander ([www.aladdinsys.com](http://www.aladdinsys.com)).

Place the folder called "CL (Classic)" anywhere you want, such as in your Applications (Mac OS 9) folder.

Launch the application by double-clicking the icon or by dragging an image onto the icon.

### Dongle installation

The dongle requires a driver to function. A driver is a small application that runs in the background.

Open the folder "Dongle installer"

Expand the archive called "hdd\_unified\_installer.sit.hqx"

Double-click the application called "aksusb Installer" and follow the instructions.

### QuickTime

QuickTime is required to run CoolLocalizer. It is developed by Apple and should already be installed on your Mac. You don't need to have the QuickTime Pro movieplayer installed, just the system extensions that probably are already installed. If QuickTime has been removed from your system you need to re-install it. Download the latest version from [www.quicktime.com](http://www.quicktime.com).

### Known issue

Certain versions of Stuffit Expander will fail to unpack the archive entirely. In this case the application icon will look like a SimpleText document called "CoolLocalizer (Classic).bin". To complete the unpacking, expand this file with Stuffit Expander.

### About dongle

A dongle is a small device required to run the full version of the application. Without the dongle installed it will run in demo/learning mode. To purchase dongles visit [www.CoolLocalizer.com](http://www.CoolLocalizer.com)

### Feedback

For support, first visit the support site at [www.coollocalizer/support.asp](http://www.coollocalizer/support.asp). Direct further support questions, suggestions, and bug reports to [support@coollocalizer.com](mailto:support@coollocalizer.com)

CoolLocalizer is ©2003 CytoLight, [www.CoolLocalizer.com](http://www.CoolLocalizer.com)