



## Windows Installation instructions

### Welcome to CoolLocalizer

CoolLocalizer© is an image analysis software for confocal microscopy.

The program was designed to be easy to use and install. Please follow the instructions below:

### Installation

Download CoolLocalizer at [www.CoolLocalizer.com](http://www.CoolLocalizer.com).

Unpack the compressed archive using Winzip or Stuffit Expander ([www.aladdinsys.com](http://www.aladdinsys.com)) by double-clicking the file.

Place the folder called "CL (Windows) folder" anywhere you want, such as on your Desktop or inside the folder "Program Files" in the C: drive. It is a good idea to place a shortcut on the desktop. To do this right-click with the mouse on the CoolLocalizer icon and drag the mouse to the desktop, select "Create shortcut here".

Launch the application by double-clicking the icon.

### Dongle installation

The dongle requires a driver to function. A driver is a small application that runs in the background.

Open the folder "Dongle installer" inside the "CL (Windows) folder".

Double-click the application called "hdd32.exe" and follow the instructions.

### QuickTime

QuickTime is required to run CoolLocalizer. It comes pre-installed on many computers and it is likely that you already have it. If CoolLocalizer fails to open images, if you get a warning by CoolLocalizer during startup, or if QuickTime has been removed from your system you need to re-install it. Download the latest version (free) from [www.quicktime.com](http://www.quicktime.com).

### About dongle

A dongle is a small device required to run the full version of the application. Without the dongle installed it will run in demo/learning mode. To purchase dongles visit [www.CoolLocalizer.com](http://www.CoolLocalizer.com)

### Feedback

For support, first visit the support site at [www.coollocalizer/support.asp](http://www.coollocalizer/support.asp). Direct further support questions, suggestions, and bug reports to [support@coollocalizer.com](mailto:support@coollocalizer.com)

CoolLocalizer is ©2003 CytoLight, [www.CoolLocalizer.com](http://www.CoolLocalizer.com)